

ETV6 Conjugated Antibody

Catalog No: #C56856



Package Size: #C56856-AF350 100ul #C56856-AF405 100ul #C56856-AF488 100ul #C56856-AF555 100ul #C56856-AF594 100ul #C56856-AF647 100ul #C56856-AF680 100ul #C56856-AF750 100ul #C56856-Biotin 100ul #C56856-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ETV6 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	ETV6 Antibody detects endogenous levels of total ETV6
Immunogen Description	A synthesized peptide derived from human ETV6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ETV6; Tel; TEL oncogene; TEL/ABL; TEL1;
Accession No.	Uniprot:P41212
Calculated MW	53kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200
IF:1:50-1:200
FC: 1:50-1:200

Product Description

Transcriptional repressor; binds to the DNA sequence 5'-CCGGAAGT-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.