

Gemin 2 Conjugated Antibody

Catalog No: #C56857

Package Size: #C56857-AF350 100ul #C56857-AF405 100ul #C56857-AF488 100ul #C56857-AF555 100ul #C56857-AF594 100ul #C56857-AF647 100ul #C56857-AF680 100ul #C56857-AF750 100ul #C56857-Biotin 100ul #C56857-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Gemin 2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Gemin 2 Antibody detects endogenous levels of total Gemin 2
Immunogen Description	A synthesized peptide derived from human Gemin 2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	gemin2; SIP 1; SIP1; SIP1 delta; SMN interacting protein 1; Survival interacting protein 1;
Accession No.	Uniprot:O14893
Calculated MW	31kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

The SMN complex plays an essential role in spliceosomal snRNP assembly in the cytoplasm and is required for pre-mRNA splicing in the nucleus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.