

CNTP2 Conjugated Antibody

Catalog No: #C56867



Package Size: #C56867-AF350 100ul #C56867-AF405 100ul #C56867-AF488 100ul #C56867-AF555 100ul #C56867-AF594 100ul #C56867-AF647 100ul #C56867-AF680 100ul #C56867-AF750 100ul #C56867-Biotin 100ul #C56867-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CNTP2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CNTP2 Antibody detects endogenous levels of total CNTP2
Immunogen Description	A synthesized peptide derived from human CNTP2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AUTS15; CDFE; CNTNAP2; CNTP2; NRXN4; PTHSL1;
Accession No.	Uniprot:Q9UHC6
Calculated MW	160kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May play a role in the formation of functional distinct domains critical for saltatory conduction of nerve impulses in myelinated nerve fibers. Seems to demarcate the juxtaparanodal region of the axo-glia junction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.