

# CBR1 Conjugated Antibody

Catalog No: #C56869

Package Size: #C56869-AF350 100ul #C56869-AF405 100ul #C56869-AF488 100ul #C56869-AF555 100ul #C56869-AF594 100ul #C56869-AF647 100ul #C56869-AF680 100ul #C56869-AF750 100ul #C56869-Biotin 100ul #C56869-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CBR1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CBR1 Antibody detects endogenous levels of total CBR1
Immunogen Description	A synthesized peptide derived from human CBR1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CBR1; CRN; SDR21C1;
Accession No.	Uniprot:P16152
Calculated MW	30kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

NADPH-dependent reductase with broad substrate specificity. Catalyzes the reduction of a wide variety of carbonyl compounds including quinones, prostaglandins, menadione, plus various xenobiotics.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.