

ERCC8 Conjugated Antibody

Catalog No: #C56871



Package Size: #C56871-AF350 100ul #C56871-AF405 100ul #C56871-AF488 100ul #C56871-AF555 100ul #C56871-AF594 100ul #C56871-AF647 100ul #C56871-AF680 100ul #C56871-AF750 100ul #C56871-Biotin 100ul #C56871-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ERCC8 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	ERCC8 Antibody detects endogenous levels of total ERCC8
Immunogen Description	A synthesized peptide derived from human ERCC8
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CKN1; CSA; ERCC8;
Accession No.	Uniprot:Q13216
Calculated MW	44kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Substrate-recognition component of the CSA complex, a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex, involved in transcription-coupled nucleotide excision repair.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.