

# MAT2A Conjugated Antibody

Catalog No: #C56884

Package Size: #C56884-AF350 100ul #C56884-AF405 100ul #C56884-AF488 100ul #C56884-AF555 100ul #C56884-AF594 100ul #C56884-AF647 100ul #C56884-AF680 100ul #C56884-AF750 100ul #C56884-Biotin 100ul #C56884-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MAT2A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	MAT2A Antibody detects endogenous levels of total MAT2A
Immunogen Description	A synthesized peptide derived from human MAT2A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AMS2; MAT 2; MAT 2A; MAT II; Mat2a; MATA2; MATII; SAMS 2; SAMS2;
Accession No.	Uniprot:P31153
Calculated MW	48kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Catalyzes the formation of S-adenosylmethionine from methionine and ATP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.