

PDE4B Conjugated Antibody

Catalog No: #C56888

Package Size: #C56888-AF350 100ul #C56888-AF405 100ul #C56888-AF488 100ul #C56888-AF555 100ul #C56888-AF594 100ul #C56888-AF647 100ul #C56888-AF680 100ul #C56888-AF750 100ul #C56888-Biotin 100ul #C56888-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PDE4B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	PDE4B Antibody detects endogenous levels of total PDE4B
Immunogen Description	A synthesized peptide derived from human PDE4B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DPDE4; PDE 4B; PDE IVB; PDE32; Pde4b; PDE4B5; PDEIVB; Phosphodiesterase 4B;
Accession No.	Uniprot:Q07343
Calculated MW	95kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes. May be involved in mediating central nervous system effects of therapeutic agents ranging from antidepressants to antiasthmatic and anti-inflammatory agents.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.