

FSH beta Conjugated Antibody

Catalog No: #C56892

Package Size: #C56892-AF350 100ul #C56892-AF405 100ul #C56892-AF488 100ul #C56892-AF555 100ul #C56892-AF594 100ul #C56892-AF647 100ul #C56892-AF680 100ul #C56892-AF750 100ul #C56892-Biotin 100ul #C56892-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	FSH beta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	FSH beta Antibody detects endogenous levels of total FSH beta
Immunogen Description	A synthesized peptide derived from human FSH beta
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FSHB; FSHbeta;
Accession No.	Uniprot:P01225
Calculated MW	20kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.