

## PCBP1 Conjugated Antibody

Catalog No: #C56893



Package Size: #C56893-AF350 100ul #C56893-AF405 100ul #C56893-AF488 100ul #C56893-AF555 100ul #C56893-AF594 100ul #C56893-AF647 100ul #C56893-AF680 100ul #C56893-AF750 100ul #C56893-Biotin 100ul #C56893-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PCBP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	PCBP1 Antibody detects endogenous levels of total PCBP1
Immunogen Description	A synthesized peptide derived from human PCBP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Alpha-CP1; hnRNP E1; HNRPE1; HNRPX; PCBP1;
Accession No.	Uniprot:Q15365
Calculated MW	40kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Single-stranded nucleic acid binding protein that binds preferentially to oligo dC.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.