

GJB1 Conjugated Antibody

Catalog No: #C56896



Package Size: #C56896-AF350 100ul #C56896-AF405 100ul #C56896-AF488 100ul #C56896-AF555 100ul #C56896-AF594 100ul #C56896-AF647 100ul #C56896-AF680 100ul #C56896-AF750 100ul #C56896-Biotin 100ul #C56896-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GJB1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	GJB1 Antibody detects endogenous levels of total GJB1
Immunogen Description	A synthesized peptide derived from human GJB1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CMTX; CMTX1; Connexin32; Cx32; GJB1;
Accession No.	Uniprot:P08034
Calculated MW	32kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.