

ATRIP Conjugated Antibody

Catalog No: #C56897



Package Size: #C56897-AF350 100ul #C56897-AF405 100ul #C56897-AF488 100ul #C56897-AF555 100ul #C56897-AF594 100ul #C56897-AF647 100ul #C56897-AF680 100ul #C56897-AF750 100ul #C56897-Biotin 100ul #C56897-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ATRIP Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	ATRIP Antibody detects endogenous levels of total ATRIP
Immunogen Description	A synthesized peptide derived from human ATRIP
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AGS1; ATIP; Atrip;
Accession No.	Uniprot:Q8WXE1
Calculated MW	86kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Required for checkpoint signaling after DNA damage. Required for ATR expression, possibly by stabilizing the protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.