

DMAP1 Conjugated Antibody

Catalog No: #C56899

Package Size: #C56899-AF350 100ul #C56899-AF405 100ul #C56899-AF488 100ul #C56899-AF555 100ul #C56899-AF594 100ul #C56899-AF647 100ul #C56899-AF680 100ul #C56899-AF750 100ul #C56899-Biotin 100ul #C56899-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	DMAP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	DMAP1 Antibody detects endogenous levels of total DMAP1
Immunogen Description	A synthesized peptide derived from human DMAP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DNA methyltransferase 1 associated protein 1; DNMAP1; DNMTAP1; EAF2; SWC4;
Accession No.	Uniprot:Q9NPF5
Calculated MW	53kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in transcription repression and activation. Its interaction with HDAC2 may provide a mechanism for histone deacetylation in heterochromatin following replication of DNA at late firing origins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.