

CC2D1A Conjugated Antibody

Catalog No: #C56901

Package Size: #C56901-AF350 100ul #C56901-AF405 100ul #C56901-AF488 100ul #C56901-AF555 100ul #C56901-AF594 100ul #C56901-AF647 100ul #C56901-AF680 100ul #C56901-AF750 100ul #C56901-Biotin 100ul #C56901-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CC2D1A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CC2D1A Antibody detects endogenous levels of total CC2D1A
Immunogen Description	A synthesized peptide derived from human CC2D1A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Coiled coil and C2 domain containing 1A; FREUD 1; Freud 1/Aki1; MRT3;
Accession No.	Uniprot:Q6P1N0
Calculated MW	130kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May play a role in the altered regulation of HTR1A associated with anxiety and major depression. It also mediates HDAC-independent repression of HTR1A promoter in neuronal cell.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.