

p66 alpha Conjugated Antibody

Catalog No: #C56903



Package Size: #C56903-AF350 100ul #C56903-AF405 100ul #C56903-AF488 100ul #C56903-AF555 100ul #C56903-AF594 100ul #C56903-AF647 100ul #C56903-AF680 100ul #C56903-AF750 100ul #C56903-Biotin 100ul #C56903-Conjugated 50ul
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	p66 alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	p66 alpha Antibody detects endogenous levels of total p66 alpha
Immunogen Description	A synthesized peptide derived from human p66 alpha
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GATAD2A; Hp66alpha; P66 alpha; p66alpha;
Accession No.	Uniprot:Q86YP4
Calculated MW	65,68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Transcriptional repressor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.