

HTF9C Conjugated Antibody

Catalog No: #C56915

Package Size: #C56915-AF350 100ul #C56915-AF405 100ul #C56915-AF488 100ul #C56915-AF555 100ul #C56915-AF594 100ul
 #C56915-AF647 100ul #C56915-AF680 100ul #C56915-AF750 100ul #C56915-Biotin 100ul #C56915-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HTF9C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	HTF9C Antibody detects endogenous levels of total HTF9C
Immunogen Description	A synthesized peptide derived from HTF9C
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HTF9C; TRMT2A;
Accession No.	Uniprot:Q8IZ69
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May be involved in nucleic acid metabolism and/or modifications.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.