

# NUDT5 Conjugated Antibody

Catalog No: #C56936



Package Size: #C56936-AF350 100ul #C56936-AF405 100ul #C56936-AF488 100ul #C56936-AF555 100ul #C56936-AF594 100ul #C56936-AF647 100ul #C56936-AF680 100ul #C56936-AF750 100ul #C56936-Biotin 100ul #C56936-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NUDT5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	NUDT5 Antibody detects endogenous levels of total NUDT5
Immunogen Description	A synthesized peptide derived from human NUDT5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hYSAH1; Nudt5; YSA1; YSA1H;
Accession No.	Uniprot:Q9UKK9
Calculated MW	30kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Hydrolyzes with similar activities ADP-ribose, ADP-mannose, and ADP-glucose. Can also hydrolyze other nucleotide sugars with low activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.