

# UBE4A Conjugated Antibody

Catalog No: #C56944

Package Size: #C56944-AF350 100ul #C56944-AF405 100ul #C56944-AF488 100ul #C56944-AF555 100ul #C56944-AF594 100ul #C56944-AF647 100ul #C56944-AF680 100ul #C56944-AF750 100ul #C56944-Biotin 100ul #C56944-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	UBE4A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	UBE4A Antibody detects endogenous levels of total UBE4A
Immunogen Description	A synthesized peptide derived from human UBE4A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	UBE4A; UBOX2; UFD2; UFD2b;
Accession No.	Uniprot:Q14139
Calculated MW	123kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.