

CHD3 Conjugated Antibody

Catalog No: #C56950

Package Size: #C56950-AF350 100ul #C56950-AF405 100ul #C56950-AF488 100ul #C56950-AF555 100ul #C56950-AF594 100ul #C56950-AF647 100ul #C56950-AF680 100ul #C56950-AF750 100ul #C56950-Biotin 100ul #C56950-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | CHD3 Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB, IF, FC |
| Species Reactivity | Human |
| Specificity | CHD3 Antibody detects endogenous levels of total CHD3 |
| Immunogen Description | A synthesized peptide derived from human CHD3 |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | Chd3; Chd7; hZFH; Mi2 ALPHA; Prp7; Prp9 1; ZFH; |
| Accession No. | Uniprot:Q12873 |
| Calculated MW | 260kDa |
| Storage | Store at 4°C in dark for 6 months |

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.