

CD58 Conjugated Antibody

Catalog No: #C56951



Package Size: #C56951-AF350 100ul #C56951-AF405 100ul #C56951-AF488 100ul #C56951-AF555 100ul #C56951-AF594 100ul
 #C56951-AF647 100ul #C56951-AF680 100ul #C56951-AF750 100ul #C56951-Biotin 100ul #C56951-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD58 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	CD58 Antibody detects endogenous levels of total CD58
Immunogen Description	A synthesized peptide derived from human CD58
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AG3; CD58; LFA3;
Accession No.	Uniprot:P19256
Calculated MW	50-75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Ligand of the T-lymphocyte CD2 glycoprotein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.