

PHD3 Conjugated Antibody

Catalog No: #C56958



Package Size: #C56958-AF350 100ul #C56958-AF405 100ul #C56958-AF488 100ul #C56958-AF555 100ul #C56958-AF594 100ul #C56958-AF647 100ul #C56958-AF680 100ul #C56958-AF750 100ul #C56958-Biotin 100ul #C56958-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PHD3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	PHD3 Antibody detects endogenous levels of total PHD3
Immunogen Description	A synthesized peptide derived from human PHD3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Egl nine homolog 3; EGLN3; Factor responsive smooth muscle protein; HIF Prolyl Hydroxylase 3; HIFP4H3; HIFPH3; P4H3; PHD3; SM20;
Accession No.	Uniprot:Q9H6Z9
Calculated MW	27kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Catalyzes the post-translational formation of 4-hydroxyproline in hypoxia-inducible factor (HIF) alpha proteins. Hydroxylates HIF-1 alpha at 'Pro-564', and HIF-2 alpha. Functions as a cellular oxygen sensor and, under normoxic conditions, targets HIF through the hydroxylation for proteasomal degradation via the von Hippel-Lindau ubiquitination complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.