

MRRF Conjugated Antibody

Catalog No: #C56965

Package Size: #C56965-AF350 100ul #C56965-AF405 100ul #C56965-AF488 100ul #C56965-AF555 100ul #C56965-AF594 100ul #C56965-AF647 100ul #C56965-AF680 100ul #C56965-AF750 100ul #C56965-Biotin 100ul #C56965-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MRRF Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	MRRF Antibody detects endogenous levels of total MRRF
Immunogen Description	A synthesized peptide derived from human MRRF
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MRFF; MTRRF; Ribosome recycling factor, mitochondrial; RRF;
Accession No.	Uniprot:Q96E11
Calculated MW	24kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.