

## GLU2B Conjugated Antibody

Catalog No: #C56966



Package Size: #C56966-AF350 100ul #C56966-AF405 100ul #C56966-AF488 100ul #C56966-AF555 100ul #C56966-AF594 100ul #C56966-AF647 100ul #C56966-AF680 100ul #C56966-AF750 100ul #C56966-Biotin 100ul #C56966-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GLU2B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GLU2B Antibody detects endogenous levels of total GLU2B
Immunogen Description	A synthesized peptide derived from human GLU2B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	G19P1; PCLD; PKCSH; PLD1; PRKCSH;
Accession No.	Uniprot:P14314
Calculated MW	80kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Regulatory subunit of glucosidase II.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.