

Fukutin Conjugated Antibody

Catalog No: #C56975



Package Size: #C56975-AF350 100ul #C56975-AF405 100ul #C56975-AF488 100ul #C56975-AF555 100ul #C56975-AF594 100ul #C56975-AF647 100ul #C56975-AF680 100ul #C56975-AF750 100ul #C56975-Biotin 100ul #C56975-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Fukutin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Fukutin Antibody detects endogenous levels of total Fukutin
Immunogen Description	A synthesized peptide derived from human Fukutin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CMD1X; FCMD; FKTN; Fukutin; LGMD2M; MDDGA4; MDDGB4; MDDGC4;
Accession No.	Uniprot:O75072
Calculated MW	51kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May be a glycosyltransferase which participates in glycosylation of alpha-dystroglycan/DAG1. May interact with and reinforce a large complex encompassing the outside and inside of muscle membranes. Could be involved in brain development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.