

GSTO1 Conjugated Antibody

Catalog No: #C56980



Package Size: #C56980-AF350 100ul #C56980-AF405 100ul #C56980-AF488 100ul #C56980-AF555 100ul #C56980-AF594 100ul
 #C56980-AF647 100ul #C56980-AF680 100ul #C56980-AF750 100ul #C56980-Biotin 100ul #C56980-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GSTO1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GSTO1 Antibody detects endogenous levels of total GSTO1
Immunogen Description	A synthesized peptide derived from human GSTO1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Gsto1; GSTTLp28; GSTX; Gtsttl; p28;
Accession No.	Uniprot:P78417
Calculated MW	28kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Exhibits glutathione-dependent thiol transferase and dehydroascorbate reductase activities. Has S-(phenacyl)glutathione reductase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.