

CPSF73 Conjugated Antibody

Catalog No: #C56984



Package Size: #C56984-AF350 100ul #C56984-AF405 100ul #C56984-AF488 100ul #C56984-AF555 100ul #C56984-AF594 100ul #C56984-AF647 100ul #C56984-AF680 100ul #C56984-AF750 100ul #C56984-Biotin 100ul #C56984-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CPSF73 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CPSF73 Antibody detects endogenous levels of total CPSF73
Immunogen Description	A synthesized peptide derived from human CPSF73
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CPSF; cpsf3; CPSF73; YSH1;
Accession No.	Uniprot:Q9UKF6
Calculated MW	77kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Component of the cleavage and polyadenylation specificity factor (CPSF) complex that play a key role in pre-mRNA 3'-end formation, recognizing the AAUAAA signal sequence and interacting with poly(A) polymerase and other factors to bring about cleavage and poly(A) addition.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.