

## Semaphorin 4D / CD100 Conjugated Antibody

Catalog No: #C57048



Package Size: #C57048-AF350 100ul #C57048-AF405 100ul #C57048-AF488 100ul #C57048-AF555 100ul #C57048-AF594 100ul #C57048-AF647 100ul #C57048-AF680 100ul #C57048-AF750 100ul #C57048-Biotin 100ul #C57048-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Semaphorin 4D / CD100 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Semaphorin 4D / CD100 Antibody detects endogenous levels of total Semaphorin 4D / CD100
Immunogen Description	A synthesized peptide derived from human Semaphorin 4D / CD100
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BB18; CD100; COLL4; Collapsin 4; GR3; MSEMA; Sema4d; Semacl2; SemaH; SEMAJ; Semaphorin 4D; Semaphorin H; Semaphorin J; Semcl2;
Accession No.	Uniprot:Q92854
Calculated MW	120,140kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF: 1:50-1:200

## Product Description

May play a functional role in the immune system, as well as in the nervous system. Induces B-cells to aggregate and improves their viability in vitro.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.