

GFPT1 Conjugated Antibody

Catalog No: #C57052

Package Size: #C57052-AF350 100ul #C57052-AF405 100ul #C57052-AF488 100ul #C57052-AF555 100ul #C57052-AF594 100ul #C57052-AF647 100ul #C57052-AF680 100ul #C57052-AF750 100ul #C57052-Biotin 100ul #C57052-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GFPT1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	GFPT1 Antibody detects endogenous levels of total GFPT1
Immunogen Description	A synthesized peptide derived from human GFPT1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CMS12; CMSTA1; GFA; GFAT; GFAT1; GFAT1m; GFPT; Gfpt1; GFPT1L; Glucosamine--fructose-6-phosphate aminotransferase [isomerizing] 1; Hexosephosphate aminotransferase 1;
Accession No.	Uniprot:Q06210
Calculated MW	79kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors for N- and O-linked glycosylation of proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.