

ICAM3 / CD50 Conjugated Antibody

Catalog No: #C57053

Package Size: #C57053-AF350 100ul #C57053-AF405 100ul #C57053-AF488 100ul #C57053-AF555 100ul #C57053-AF594 100ul #C57053-AF647 100ul #C57053-AF680 100ul #C57053-AF750 100ul #C57053-Biotin 100ul #C57053-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ICAM3 / CD50 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	ICAM3 / CD50 Antibody detects endogenous levels of total ICAM3 / CD50
Immunogen Description	A synthesized peptide derived from human ICAM3 / CD50
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD50; CDw50; ICAM 3; ICAM3; ICAMR; Intercellular adhesion molecule 3;
Accession No.	Uniprot:P32942
Calculated MW	100-160kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). ICAM3 is also a ligand for integrin alpha-D/beta-2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.