

EGFR (Phospho-Tyr845) Conjugated Antibody

Catalog No: #C14294

Package Size: #C14294-AF350 100ul #C14294-AF405 100ul #C14294-AF488 100ul #C14294-AF555 100ul #C14294-AF594 100ul #C14294-AF647 100ul #C14294-AF680 100ul #C14294-AF750 100ul #C14294-Biotin 100ul #C14294-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EGFR (Phospho-Tyr845) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-EGFR (Y845) Antibody detects endogenous levels of total Phospho-EGFR (Y845)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EGFR; ERBB; ERBB1; Errp; HER1; mENA; NISBD2; PIG61; SA7; wa2; Wa5;
Accession No.	P00533
Calculated MW	175kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

Receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.