

## AKT1/2/3 Conjugated Antibody

Catalog No: #C57129



Package Size: #C57129-AF350 100ul #C57129-AF405 100ul #C57129-AF488 100ul #C57129-AF555 100ul #C57129-AF594 100ul #C57129-AF647 100ul #C57129-AF680 100ul #C57129-AF750 100ul #C57129-Biotin 100ul #C57129-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AKT1/2/3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	AKT1 + AKT2 + AKT3 Antibody detects endogenous levels of total AKT1 + AKT2 + AKT3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PKB;RAC;PRKBA;PKB-ALPHA;RAC-ALPHA; AKT1 ; AKT2; AKT3; Akt1/2/3; pan Akt;
Accession No.	P31749/P31751/Q9Y243
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Background

A synthesized peptide derived from human AKT1 + AKT2 + AKT3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.