

HN1 Conjugated Antibody

Catalog No: #C57132

Package Size: #C57132-AF350 100ul #C57132-AF405 100ul #C57132-AF488 100ul #C57132-AF555 100ul #C57132-AF594 100ul #C57132-AF647 100ul #C57132-AF680 100ul #C57132-AF750 100ul #C57132-Biotin 100ul #C57132-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	HN1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	HN1 Antibody detects endogenous levels of total HN1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ARM2; HN1; HN1A;
Accession No.	Q9UK76
Calculated MW	21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Background

Modulates negatively AKT-mediated GSK3B signaling (PubMed:21323578, PubMed:22155408).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.