

CD223 Conjugated Antibody

Catalog No: #C57166

Package Size: #C57166-AF350 100ul #C57166-AF405 100ul #C57166-AF488 100ul #C57166-AF555 100ul #C57166-AF594 100ul #C57166-AF647 100ul #C57166-AF680 100ul #C57166-AF750 100ul #C57166-Biotin 100ul #C57166-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD223 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	IF
Species Reactivity	Human
Specificity	CD223 Antibody detects endogenous levels of total CD223
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD223; Lag3; Lymphocyte activating 3; Protein FDC;
Accession No.	P18627
Calculated MW	60-80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Background

Involved in lymphocyte activation. Binds to HLA class-II antigens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.