

Human FGF19/ Fibroblast growth factor 19 ELISA Kit

Catalog No: #EK21106



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FGF19/ Fibroblast growth factor 19 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:7.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.