

# PE mouse IgG2a isotype Control

Catalog No: #IC002aPE



Package Size: #IC002aPE-1 25tests #IC002aPE-2 100tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE mouse IgG2a isotype Control
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4AC44
Isotype	Mouse IgG2a
Purification	Purified by Protein A.
Applications	FC IF
Concentration	10 <sup>3</sup> OI/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA, Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.