

Human Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit

Catalog No: #EK20505

Package Size: #EK20505-1 48T #EK20505-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range: 1.25-80 ng/mL

Sensitivity: 0.48 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.