

Goat Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK20521



Package Size: #EK20521-1 48T #EK20521-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goat Immunoglobulin G (IgG) ELISA Kit
Brief Description	ELISA Kit
Host Species	Goat
Applications	ELISA

Application Details

Detect Range: 9.37-600 µg/mL

Sensitivity: 3.5 µg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.