

# ZO-1 Conjugated Antibody

Catalog No: #C29274

Package Size: #C29274-AF350 100ul #C29274-AF405 100ul #C29274-AF488 100ul #C29274-AF555 100ul #C29274-AF594 100ul  
 #C29274-AF647 100ul #C29274-AF680 100ul #C29274-AF750 100ul #C29274-Biotin 100ul #C29274-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZO-1 Conjugated Antibody
Host Species	Rabbit
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZO-1 (NP_003248.3).
Other Names	TJP1;ZO-1;ZO1
Calculated MW	195kDa/250kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF: 1:50-1:200

## Background

This gene encodes a protein located on a cytoplasmic membrane surface of intercellular tight junctions. The encoded protein may be involved in signal transduction at cell-cell junctions. Alternative splicing of this gene results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.