

# Human IgG Conjugated Antibody

Catalog No: #CIG701



Package Size: #CIG701-AF350 100ul #CIG701-AF405 100ul #CIG701-AF488 100ul #CIG701-AF555 100ul #CIG701-AF594 100ul #CIG701-AF680 100ul #CIG701-AF750 100ul #CIG701-Biotin 100ul #CIG701-F647 100ul  
Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human IgG Conjugated Antibody
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,95%
Concentration	1 mg/ml
Storage	Store at 4°C

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
AF405 conjugated: most applications: 1: 50 - 1: 250
AF488 conjugated: most applications: 1: 50 - 1: 250
AF555 conjugated: most applications: 1: 50 - 1: 250
AF594 conjugated: most applications: 1: 50 - 1: 250
AF647 conjugated: most applications: 1: 50 - 1: 250
AF680 conjugated: most applications: 1: 50 - 1: 250
AF750 conjugated: most applications: 1: 50 - 1: 250
Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.