

CD69 Monoclonal Antibody

Catalog No: #48135

Package Size: #48135-1 50ul #48135-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD69 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification
Applications	WB;IHC;IF;FCM;ELISA
Species Reactivity	Human
Specificity	CD69 Monoclonal Antibody detects endogenous levels of CD69 protein.
Immunogen Description	Purified recombinant fragment of human CD69 expressed in E. Coli.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Storage	Store at -20°C/1 year

Application Details

WB: 1/500 - 1/2000.


IHC: 1/200 - 1/1000.

Flow cytometry: 1/200 - 1/400.


ELISA: 1/10000.

IF 1:50-200

Images



Western Blot analysis using CD69 Monoclonal Antibody against, Jurkat (1), L1210 (2) and TPH-1 (3) cell lysate.



Immunohistochemistry analysis of paraffin-embedded human Tonsil tissues with AEC staining using CD69 Monoclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.