

Monkey N-terminal pro-brain natriuretic peptide (NT-proBNP) ELISA Kit

Catalog No: #EK20523



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20523-1 48T #EK20523-2 96T

Description

Product Name	Monkey N-terminal pro-brain natriuretic peptide (NT-proBNP) ELISA Kit
Brief Description	Monkey N-terminal pro-brain natriuretic peptide (NT-proBNP) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:1.56-100ng/mL
Sensitivity:0.62ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.