

β-Actin Rabbit Conjugated Antibody

Catalog No: #C52901

Package Size: #C52901-AF350 100ul #C52901-AF405 100ul #C52901-AF488 100ul #C52901-AF555 100ul #C52901-AF594 100ul #C52901-AF647 100ul #C52901-AF680 100ul #C52901-AF750 100ul #C52901-Biotin 100ul #C52901-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	β-Actin Rabbit Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	beta Actin Antibody detects endogenous levels of total beta-Actin protein.
Immunogen Description	A synthesized peptide derived from human beta Actin
Other Names	ACTIN; ACTB; BETA-ACTIN; beta actin; actin beta;
Calculated MW	42kDa
Formulation	Liquid
Storage	Store at 4Λ C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.