

β -Tubulin Rabbit Conjugated Antibody

Catalog No: #C52905



Package Size: #C52905-AF350 100ul #C52905-AF405 100ul #C52905-AF488 100ul #C52905-AF555 100ul #C52905-AF594 100ul #C52905-AF647 100ul #C52905-AF680 100ul #C52905-AF750 100ul #C52905-Biotin 100ul #C52905-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	β -Tubulin Rabbit Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat Cow Dog Chicken Monkey Zebrafish Xenopus tropicalis Drosophila melanogaster
Specificity	beta Tubulin Antibody detects endogenous levels of total beta Tubulin
Immunogen Description	A synthesized peptide derived from human beta Tubulin
Other Names	Beta 4 tubulin; Beta 5 tubulin; Beta1 tubulin; Tubulin beta chain; Tubulin beta class I;
Calculated MW	50kDa
Formulation	Liquid
Storage	Store at 4 Λ C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.