

Recombinant Human Chromobox protein homolog 5(CBX5)

Catalog No: #AP89655

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89655-1 20ug #AP89655-2 100ug #AP89655-3 1mg

Description

Product Name	Recombinant Human Chromobox protein homolog 5(CBX5)
Brief Description	Recombinant protein
Host Species	E.coli
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Protein Length:Full Length Expression Region: 1-191aa
Other Names	Antigen p25; CBX5; CBX5_HUMAN; CG8409; Chromobox 5; Chromobox homolog 5 (HP1 alpha homolog; Drosophila); Chromobox homolog 5; Chromobox protein homolog 5; Epididymis luminal protein 25; HEL25; Heterochromatin protein 1 alpha; Heterochromatin protein 1; Heterochromatin protein 1 homolog alpha; HP1 alpha; HP1 alpha homolog; HP1; HP1A; HP1Hs alpha; Su(var)205
Target Species	Homo sapiens (Human)
Calculated MW	29.7 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	MGKKTkRTADSSSEDEEEYVVEKVLDRRVKGGQVEYLLKWKGFSEEHNtWEPEKNLDCPELISEFMKKYKK MKEGENNkPREKSESnkRKSnfSNSADDIKSKKKREQSNDIARGFERGLEPEKIIGATDSCGDLMLFMKWKD TDEADLVLAKEANVKCPQIVIAFYEERLTWHAYPEDAENKEKETAKS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.