

Recombinant Human Delta-type opioid receptor(OPRD1)

Catalog No: #AP89656

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89656-1 20ug

Description

Product Name	Recombinant Human Delta-type opioid receptor(OPRD1)
Brief Description	Recombinant Human Delta-type opioid receptor(OPRD1)
Host Species	E.coli
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Protein Length:Full Length Expression Region: 1-372aa
Target Name	His-tagged
Target Species	Homo sapiens (Human)
Calculated MW	41.9 kDa
Formulation	Protein Length:Full Length Expression Region: 1-372aa MEPAPSAGAELQPPLFANASDAYPSACPSAGANASGPPGARSASSLALAIITALYSAVCAVGLLGNVLMVFGI VRYTKMKTATNIYIFNLALADALATSTLPFQSAKYLMEWPFGELLCKAVLSIDYYNMFTSIFTLTMMMSVDRYIAV CHPVKALDFRTPAKAKLINICIWVLASGVGVPIMVMAVTRPRDGAVVCMQLQFPSPSWYWDVTVKICVFLFAFVV PILIIIVCYGLMLLRSLRSVRLLSGSKEKDRSLRRITRMVLVVVGAFFVVCWAPIHIFVIVWTLVDIDRRDPLVVAALH LCIALGYANSSLNPNVLYAFLDENFKRCFRQLCRKPCGRPDPSFSRAREATARERV TACTPSDGPGGGAAA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.