

Mouse Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit

Catalog No: #EK20526

Package Size: #EK20526-1 48T #EK20526-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range: 1.25-80 mU/mL

Sensitivity: 0.49 mU/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.