

# Human Deoxycorticosterone (DOC) ELISA Kit

Catalog No: #EK20541



Package Size: #EK20541-1 48T #EK20541-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Deoxycorticosterone (DOC) ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range: 3.9 pg/ml - 250 pg/ml

Sensitivity: 1.5 pg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Detection Method: competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.