

# Human 11-Deoxycortisol ELISA Kit

Catalog No: #EK20527



Package Size: #EK20527-1 48T #EK20527-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human 11-Deoxycortisol ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:1.25pmol/mL-20 pmol/mL

Sensitivity:0.625pmol/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.