

Versican Conjugated Antibody

Catalog No: #C48600

Package Size: #C48600-AF350 100ul #C48600-AF405 100ul #C48600-AF488 100ul #C48600-AF555 100ul #C48600-AF594 100ul #C48600-AF647 100ul #C48600-AF680 100ul #C48600-Biotin 100ul #C48600-Conjugated 50ul #C48600-AF750 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Versican Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB, IF
Species Reactivity	Human; Mouse;Rat
Specificity	Versican Polyclonal Antibody detects endogenous levels of Versican protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human VCAN.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	The antiserum was produced against synthesized peptide derived from human VCAN.
Calculated MW	373kD
Concentration	1mg/mL
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.