

SSCηO ε ²ζA²(20x,pH7.0)

Catalog No: #SA7803



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name SSCηO ε ²ζA²(20x,pH7.0)

Application Details

SSCηO ε ²ζA²(20x,pH7.0) is a highly sensitive and specific antibody used for the detection of RNA and DNA. It is suitable for Northern blotting, Southern blotting, and in situ hybridization. The antibody is conjugated with a fluorophore, allowing for the visualization of the target nucleic acid. It is used in a variety of applications, including the detection of specific genes and the study of gene expression patterns. The antibody is stable and easy to use, making it a valuable tool for researchers in molecular biology and genetics.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.