

## NRIP1 Rabbit pAb

Catalog No: #55653

Package Size: #55653-1 50ul #55653-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NRIP1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant fusion protein of human NRIP1(NP_003480.2).
Other Names	NRIP1;RIP140
Calculated MW	126kDa
SDS-PAGE MW	150KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB, 1:500 - 1:2000

## Images



Western blot analysis of extracts of Mouse brain, using NRIP1 Polyclonal Antibody at 1:1000 dilution.

## Background

Nuclear receptor interacting protein 1 (NRIP1) is a nuclear protein that specifically interacts with the hormone-dependent activation domain AF2 of nuclear receptors. Also known as RIP140, this protein modulates transcriptional activity of the estrogen receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.